

SIZE CHART

Cervical Aids & Splints

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Spondylosis Collar (Philadelphia) Ethafoam	Measure distance from lower jaw to clavicle bone while in extension position.			10-11	11-12			
Spondylosis Collar (Philadelphia) Plastazote			9-10	10-11	11-12			
Soft Cervical Collar			9-10	10-11	11-12	12-14		
Soft Cervical Collar with Eyelet			9-10	10-11	11-12	12-14		
Emergency Collar			9-10	10-11	11-12	12-14		
Postural Collar (ESL)			9-10	10-11	11-12	12-14		
Hard Collar (Adjustable)			9-10	10-11	11-12			
Post Procedural Collar			9-10	10-11	11-12			
Cervical Collar (Adjustable)		7-8	9-14					
Cervical Orthosis	Measure Neck Circumference.	25-33	33-40					

Back & Thoracic Supports

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Lumbo Sacral Support	Measure circumference of the waist below the navel.		66-76	79-89	91-102	104-114	117-127	
Back Support Executive			66-76	79-89	91-102	104-114	117-127	
Front Line Back Support			66-76	79-89	91-102	104-114	117-127	
Rib Belt	Measure around chest at the base the sternum.		71-86	89-104	107-122	124-142		
Hyper Extension Brace (ASH Brace)	Measure height of the patient		up to 152		above 152			
Dorso Lumbar Spinal Orthosis (Taylor's Brace)	Measure circumference of the waist below the navel.		66-76	79-89	91-102	104-114	117-127	
Sternal Binder	Measure around chest at the base of the sternum.		71-86	86-102	102-117	117-132		
Abdominal Binder (Available in 8", 9" & 10" wide length.)	Measure circumference of the waist below the navel.		66-76	79-89	91-102	104-114	117-127	
Scrotal Support			66-76	79-89	91-102	104-114	117-127	
Re-Enforced L.S Belt			66-76	79-89	91-102	104-114	117-127	
Lumbo Sacral Pull Back Orthosis			70-100		100-130			
Hernia Belt			66-76	79-89	91-102	104-114	117-127	
Maternity Belt			86-97	97-107	107-117	117-127	127-137	
Postural Brace			66-76	79-89	91-102	104-114	117-127	
Comfort Back Support			66-76	79-89	91-102	104-114	117-127	
Lumbar Support								66-127
Spinal Brace Adjustable								66-127

Shoulder, Arm & Clavicle Supports

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Shoulder Immobilizer with Shoulder Support	Measure around chest at the base of the sternum.		71-86	89-104	107-122	124-142		
Arm Sling Pouch	Measure forearm from elbow to knuckles		18-23	25-30	33-38	41-46		
Arm Sling Strap								18-46
Clavicle Strap	Measure around chest at the base of the sternum.		46-51	53-64	66-76			
Clavicle Support			70-85	85-100	100-115	115-132		
Humerus Support	Measure circumference at armpit		25-31	33-38	41-46			
Universal Shoulder Immobilizer	Measure forearm from elbow to knuckles.		18-23	25-30	33-38	41-46		
Shoulder Abduction Pouch			18-23	25-30	33-38	41-46		
Cuf & Collar Sling								18-46
Arm Sling Pouch (Perforated)			18-23	25-30	33-38	41-46		

Wrist, Elbow & Arm Splints

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Wrist & Forearm Splint	Measure circumference above wrist joint		15-18	20-23	25-28			
Wrist Cock – Up Splint			15-18	20-23	25-28			
Elbow Brace	Measure upper arm around the middle of biceps		23-28	30-36	38-48			
ROM Elbow Brace			23-28	30-36	38-48			
Tennis Elbow Support	Measure circumference just below the elbow joint, with the elbow bent		18-23	25-30	33-38			
Wrist Wrap	Measure circumference around the wrist joint							15-28
Wrist Wrap & Thumb Support								15-28
Thumb Spica Splint								15-28
Day Splint								15-28
Night Splint								15-28
Carpal Tunnel Splint								15-28
Dynamic Cockup Splint (Left/Right)	Measure forearm from elbow to knuckles							18-46
Hi-Tech Elbow Support	Measure circumference around elbow joint		18-23	25-30	33-38			
Hi-Tech Elbow Brace with Gel Pad			18-23	25-30	33-38			
Hi-Tech Wrist Band	Measure circumference around the wrist joint		15-18	20-23	25-28			
Hi-Tech Wrist & Forearm			15-18	20-23	25-28			
Elbow Support (Multi Directional Stretch)	Measure circumference around elbow joint		18-23	25-30	33-38			
Wrist Support Lace-Pull	Measure circumference above wrist joint		15-20		25-30			
Wrist Hand Resing Splint - Finger	Measure width of metacarpal joints (1-4)		8-10	10-13	13-15			
Wrist Hand Resing Splint – Cone			8-10	10-13	13-15			

Knee & Leg Supports

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Elastic Knee Support with hinges	Measure circumference 3" above patella		33-38	41-46	48-53	56-61		
Knee Immobilizer 12"			33-38	41-46	48-53	56-61		
Knee Immobilizer 14"			33-38	41-46	48-53	56-61		
Knee Immobilizer 19"			33-38	41-46	48-53	56-61		
Knee Immobilizer 22"			33-38	41-46	48-53	56-61		
Knee Immobilizer Universal								33-61
ROM Knee Brace			33-38	41-46	48-53	56-61		
Knee Cap (Multi Directional Stretch)			33-38	41-46	48-53	56-61	63-68	
Femur Support	Measure circumference at mid - thigh.		46-51	53-58	61-69			
Tibial Support	Measure circumference at calf.		25-31	33-38	41-46			
Knee Cap Open Patella with Gel Ring	Measure circumference 3" above patella		33-38	41-46	48-53	56-61	63-68	
Hi-Tech Elastic Knee Support			33-38	41-46	48-53	56-61	63-68	
Hi-Tech Knee Support with Hinges & Gel Pad			33-38	41-46	48-53	56-61	63-68	
Knee Support Octopus Type with Hinge			33-38	41-46	48-53	56-61	63-68	
ROM Knee Brace Universal								33-61

Ankle & Foot Splints

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Ankle Foot Orthosis (Left / Right)	Measure foot size	10-15	18-20	23-25	28-30			
Foot Drop Splint (Left/ Right)	Measure length from heel towards calf.	10-15	18-23	25-30	33-38	41-46		
Cast Shoe	Measure foot size	10-15	18-20	23-25	28-30			
De Rotation Boot	Measure foot size	10-15						18-30
Ankle Brace	Measure circumference just above ankle joint		15-18	20-23	25-28			
Ankle Support								15-28
Anklet			15-18	20-23	25-28			
Ankle Wrap			15-18	20-23	25-28			
Anklet with Ankle Wrap			15-18	20-23	25-28			
Hi-Tech Ankle Support			15-18	20-23	25-28			
Hi-Tech Ankle Support with Gel Pad			15-18	20-23	25-28			
Air Walker Pneumatic (Short / Long)	Measure foot size		18-20	23-25	28-30			

Finger Splints

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Cot Splint	Measure length of finger		Below 5	5-7.6	7.6-10			
Frog Splint			Below 5	5-7.6	7.6-10			
Baseball Splint			Below 5	5-7.6	7.6-10			
Stax Mallet			Below 5	5-7.6	7.6-10			
Protector Splint			Below 5	5-7.6	7.6-10			
Spoon Splint			Below 5	5-7.6	7.6-10			
Extension Splint			Below 5	5-7.6	7.6-10			
Flexion Splint			Below 5	5-7.6	7.6-10			
Finger Exerciser								5-10

Neoprene Products

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Knee Support with Hinges	Measure circumference 3" above patella.		33-38	41-46	48-53	56-61	63-68	
Knee cap Pull on (Open Patella)			33-38	41-46	48-53	56-61	63-68	
Ankle Wrap	Measure circumference just above ankle joint.		15-18	20-23	25-28			
Wrist Wrap	Measure circumference around the wrist joint.							15-28
Elbow Support	Measure circumference around elbow joint.		18-23	25-30	33-38	40-45		
Wrist & Thumb Support	Measure circumference around the wrist joint.							15-28
Shoulder Support	Measure around chest at the base of the sternum.		71-86	89-104	107-122	124-142		
Wrist & Forearm brace	Measure circumference around the wrist joint		15-18	20-23	25-28			
Lumbo Sacral Belt with Gel Pad	Measure circumference of the waist below the navel.		66-76	79-89	91-102	104-114	117-127	
Knee Support Open Patella with Gel Ring	Measure circumference 3" above patella.		33-38	41-46	48-53	56-61	63-68	
Knee Support Open Patella (Octopus Type)			33-38	41-46	48-53	56-61	63-68	
Ankle Support	Measure circumference just above ankle joint		15-18	20-23	25-28			
OA Knee Support	Measure circumference 3" above patella.		33-38	41-46	48-53	56-61	63-68	

Traction Kits

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Leg Traction Kit	For Traction Pad, Measure circumference at calf.		25-31	33-38	41-46			
Skin Traction Kit Elastic			25-31	33-38	41-46			
Cervical Traction with Head Halter Sitting	Measure distance from lower jaw to clavicle bone while in extension position.		9-10	10-11	11-12			
Cervical Traction Head Halter			9-10	10-11	11-12			
Weight Cuf (1/2kg, 1kg, 2kg, 3kg)	Universal size, fits all Wrist & Ankles.							15-28
Pelvic Traction Kit	Measure circumference of the waist below the navel.		66-76	79-89	91-102	104-114	117-127	
Skin Traction Pad Crepe	Measure distance from knee to foot.	Below 30						Above 30

Silicon Products

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Heel Cushion	Measure foot size		18-20	23-25	28-30			
Insole			18-20	23-25	28-30			
Arch Pad			18-20	23-25	28-30			
Metatarsal Pad	Measure across the width at the ball of foot region		3-4	5-6	8-9			
Gel Ball	Available in Ortho (Hard) & Neuro (Soft) with Diameter of 4cm, 5cm & 6cm							
Toe Seperator	Measure the length of the inner side of the toes to be spread.		1.7	2.1	2.5			


Resuscitation Splints

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Short Arm Splint	Measure length of injured arm part to be splinted.		30.48 x 7.62					
Full Arm Splint					53.34 x 7.62			
Short Leg Splint				45.72 x 10.2				
Full Leg Splint						71.12 x 10.2		
Stockinette	Widths available: 5 cm, 7.5 cm, 10 cm, and 12.5 cm, all offered in 10 - meter lengths.							
Bunion / Hallux Valgus Splint	Measure foot size							18-30
Crepe Bandage 10cm & 15cm	Widths available: 6 cm, 8 cm, 10 cm, and 15 cm, each with a stretched length of 4 meters.							
Heating Belt	Available size: 27.9 cm x 20.3 cm							
Exercise Ball	Available in Neuro (Soft) with Diameter of 6.5cm							
Air Splints	Available in Set of 6 Splints							

Hot & Cold Packs

PRODUCT	MEASUREMENT (in CM)	Child	Small	Medium	Large	X-Large	XX-Large	Universal
Hot & Cold Compression - Wrist Wrap	Measure circumference around the wrist joint							15-28
Hot & Cold Compression - Shoulder Wrap	Measure around chest at the base of the sternum.							71-142
Hot & Cold Compression - Knee Wrap	Measure circumference 3" above patella.							33-68
Hot & Cold Compression - Ankle Wrap	Measure circumference just above ankle joint							15-28
Hot & Cold Pack	Available size: 30.4 cm x 17.7 cm & 25.4 cm x 10.16 cm							

Compression Garments

PRODUCT	MEASUREMENT (in CM)				
Graduated compression levels 15-20mmHg/ 20-30mmHg/ 30-40mmHg	Size		b	c	G
Medical Compression Stocking Knee High Class 2	<div>Refer image for stockings measurement.</div> 	Small	20 – 23	30 – 36	44 – 55
Medical Compression Stocking Thigh High Class 2		Medium	23 – 26	34 – 40	50 – 61
		Large	26 - 29	38 – 44	56 – 68
Anti Embolism Stockings Class 1 Knee High		X- Large	29 – 32	42 – 48	62 – 70
		XX- Large	32 – 35	46 – 52	68 - 82
Anti Embolism Stockings Class 1 Thigh High		Length			
			Short	Regular	Long
	/D	< 36cm	36 – 40 cm	>40 cm	
	/G	< 70cm	70 – 80 cm	>80 cm	

Body Shapers

PRODUCT	MEASUREMENT (in CM)	Small	Medium	Large	X-Large	XX-Large	XXX-Large
Body Shaper – Slip On	Measure circumference of the waist below the navel.	65 – 70	75 – 80	85 – 90	95 – 100	105 – 110	115 – 120
Body Shaper- Short Type		65 – 70	75 – 80	85 – 90	95 – 100	105 – 110	115 – 120
Body Shaper – Panty Style		65 – 70	75 – 80	85 – 90	95 – 100	105 – 110	115 – 120
Body Shaper – Fit Corset		65 – 70	75 – 80	85 – 90	95 – 100	105 – 110	115 – 120
Body Shaper – Hip Corset		65 – 70	75 – 80	85 – 90	95 – 100	105 – 110	115 – 120

Kids Range

PRODUCT	MEASUREMENT (in CM)	Small	Medium	Large
Soft Collar	Measure circumference of the neck.	17.5 - 20	20 - 22.5	22.5 - 25
Arm Sling Pouch	Measure forearm from elbow to knuckles	15 - 18	18 - 21	21 - 24
Wrist Forearm Support	Measure circumference around the wrist joint	7 - 9.5	9.5 - 11.4	11.4 - 14
Clavicle Support	Measure circumference of the chest	40 - 45	45 - 50	50 - 55
Knee Immobilizer	Measure circumference 4" above patella.	25 - 30	30 - 35	35 - 40
Cast Shoe	Measure foot size.	8 - 15		